

Specification

• Double Pigtail Drainage Catheter

| 모델명 NO. | O.D (Fr) | Pigtail (cm) | between Pigtail and Pigtail Length (cm) | Total length (cm) |
|-------------|-------------|-----------------|--|----------------------|
| EQDP071-1-1 | 7 | 1cm-1cm | 1 | 30 |
| EQDP071-2-1 | 7 | 1cm-1cm | 2 | 30 |
| EQDP072-6-1 | 7 | 2cm-1cm | 6 | 30 |
| EQDP081-1-1 | 8.5 | 1cm-1cm | 1 | 34 |
| EQDP081-2-1 | 8.5 | 1cm-1cm | 2 | 34 |
| EQDP082-1-1 | 8.5 | 2cm-1cm | 1 | 34 |
| EQDP101-1-1 | 10.2 | 1cm-1cm | 1 | 34 |
| EQDP101-2-1 | 10.2 | 1cm-1cm | 2 | 34 |
| EQDP102-1-1 | 10.2 | 2cm-1cm | 1 | 34 |

• Straight : All Purpose Drainage Catheter

| 모델명 NO. | O.D (Fr) | Total length (cm) |
|------------|-------------|----------------------|
| EQDA0600 | 6 | 34 |
| EQDA0700 | 7 | 34 |
| EQDA0800 | 8.5 | 34 |
| EQDA1000 | 10.2 | 34 |
| EQDA1200 | 12 | 34 |
| EQDA1400 | 14 | 34 |
| EQDA1600 | 16 | 34 |
| EQCDA0700 | 7 | 30 |
| EQCDA0800 | 8.5 | 30 |
| EQCDA1000 | 10.2 | 30 |
| EQCDA1200 | 12 | 30 |
| EQCDA1400 | 14 | 30 |

• Straight : Nephrostomy Drainage Catheter

| 모델명 NO. | O.D (Fr) | Total length (cm) |
|------------|-------------|----------------------|
| EQDN0500 | 5 | 34 |
| EQDN0600 | 6 | 34 |
| EQDN0700 | 7 | 34 |
| EQDN0800 | 8.5 | 34 |
| EQDN1000 | 10.2 | 34 |
| EQDN1200 | 12 | 34 |
| EQDN1400 | 14 | 34 |
| EQDN1600 | 16 | 34 |
| EQCDN0800 | 8.5 | 30 |
| EQCDN1000 | 10.2 | 30 |
| EQCDN1200 | 12 | 30 |
| EQCDN1400 | 14 | 30 |
| EQCDN1600 | 16 | 30 |

Distributor



Intervention Medical Co.,Ltd.

42, Seongnae-ro 6-gil, Gangdong-gu (Seongnae-dong, jungang Bldg, Rm202) Seoul, 05398 Republic of Korea
Tel. +82-2-747-4594 Fax. +82-2-747-4597
blueleopard@intervention-m.com

Manufacturer



HEPATO Inc.

42, Seongnae-ro 6-gil, Gangdong-gu (Seongnae-dong, jungang Bldg, Rm302) Seoul, 05398 Republic of Korea
Tel. +82-2-446-3586 Fax. +82-2-423-0749



VASFLEX Co., Ltd.

#603, #604, #605, #606, #703, 76, Dongnae-ro, Dong-gu, Daegu, Republic of Korea
Tel. +82-53-964-4249 Fax. +82-53-964-4245




EQ

DRAINAGE CATHETER




HEPATO
MEDICAL DEVICE


EQ DRAINAGE CATHETER



Safety Push Lock System



Kink Proof



Small Size Pigtail Bigger Size Hole



A double pigtail available with various hole-size options

Specification

• All Purpose Drainage Catheter

| 모델명 NO. | O.D (Fr) | Pigtail (cm) | Total length (cm) |
|------------------------|-------------|-----------------|----------------------|
| EQDA0510 | 5 | 1 | 30 |
| EQDA0520 | 5 | 2 | 30 |
| EQDA0610 | 6 | 1 | 30 |
| EQDA0620 | 6 | 2 | 30 |
| EQDA070 | 7 | 0.5 | 30 |
| EQDA0710/ EQCDA0710 | 7 | 1 | 30 |
| EQDA0720/ EQCDA0720 | 7 | 2 | 30 |
| EQDA0805 | 8.5 | 0.5 | 34 |
| EQDA0810 | 8.5 | 1 | 34 |
| EQDA0820 | 8.5 | 2 | 34 |
| EQDA1010 | 10.2 | 1 | 34 |
| EQDA1020 | 10.2 | 2 | 34 |
| EQDA1220 | 12 | 2 | 34 |
| EQDA1420 | 14 | 2 | 34 |
| EQDA1620 | 16 | 2 | 34 |
| EQCDA0805 | 8.5 | 0.5 | 30 |
| EQCDA0810 | 8.5 | 1 | 30 |
| EQCDA0820 | 8.5 | 2 | 30 |
| EQCDA1010 | 10.2 | 1 | 30 |
| EQCDA1020 | 10.2 | 2 | 30 |
| EQCDA1220 | 12 | 2 | 30 |
| EQCDA1420 | 14 | 2 | 30 |

• Biliary Drainage Catheter

| 모델명 NO. | O.D (Fr) | Pigtail (cm) | Total length (cm) |
|------------|-------------|-----------------|----------------------|
| EQDB0720 | 7 | 2 | 40 |
| EQDB0810 | 8.5 | 1 | 40 |
| EQDB0820 | 8.5 | 2 | 40 |
| EQDB1010 | 10.2 | 1 | 40 |
| EQDB1020 | 10.2 | 2 | 40 |
| EQDB1220 | 12 | 2 | 40 |
| EQDB1420 | 14 | 2 | 40 |
| EQDB1620 | 16 | 2 | 40 |

• Nephrostomy Drainage Catheter

| 모델명 NO. | O.D (Fr) | Pigtail (cm) | Total length (cm) |
|------------------------|-------------|-----------------|----------------------|
| EQDN0510 | 5 | 1 | 30 |
| EQDN0520 | 5 | 2 | 30 |
| EQDN0610 | 6 | 1 | 30 |
| EQDN0620 | 6 | 2 | 30 |
| EQDN0705/ EQCDN0705 | 7 | 0.5 | 30 |
| EQDN0710/ EQCDN0710 | 7 | 1 | 30 |
| EQDN0720/ EQCDN0720 | 7 | 2 | 30 |
| EQDN0805 | 8.5 | 0.5 | 34 |
| EQDN0810 | 8.5 | 1 | 34 |
| EQDN0820 | 8.5 | 2 | 34 |
| EQDN1010 | 10.2 | 1 | 34 |
| EQDN1020 | 10.2 | 2 | 34 |
| EQDN1220 | 12 | 2 | 34 |
| EQDN1420 | 14 | 2 | 34 |
| EQDN1620 | 16 | 2 | 34 |
| EQCDN0810 | 8.5 | 1 | 30 |
| EQCDN0820 | 8.5 | 2 | 30 |
| EQCDN1010 | 10.2 | 1 | 30 |
| EQCDN1020 | 10.2 | 2 | 30 |
| EQCDN1220 | 12 | 2 | 30 |
| EQCDN1420 | 14 | 2 | 30 |

| 모델명 NO. | O.D (Fr) | Pigtail (cm) | Total length (cm) |
|------------|-------------|-----------------|----------------------|
| EQCDB0710 | 7 | 1 | 47 |
| EQCDB0720 | 7 | 2 | 47 |
| EQCDB0810 | 8.5 | 1 | 47 |
| EQCDB0820 | 8.5 | 2 | 47 |
| EQCDB1010 | 10.2 | 1 | 47 |
| EQCDB1020 | 10.2 | 2 | 47 |
| EQCDB1220 | 12 | 2 | 47 |
| EQCDB1420 | 14 | 2 | 47 |
| EQCDB1620 | 16 | 2 | 47 |