



*Tokai Occlusion Balloon Catheter*

# **RESCUE BALLOON**

*Occlusion Catheter*


# RESCUE BALLOON

Designed Exclusively For IABO-REBOA

## 7Fr Sheath Compatible

New IABO catheter with the smallest balloon profile provides less invasive and much safety treatment in emergency cases.

## High Compliant Balloon



$\phi 16\text{mm}$

$\phi 40\text{mm}$

### SPECIFICATION

Product Code	Nominal Inflated Diameter	Shaft Diameter	Usable Length	Balloon Length	Minimum Sheath	Maximum Guidewire
RB-167080-E	16mm	7Fr	80cm	60mm	7Fr	0.025inch



**TOKAI MEDICAL PRODUCTS**

1485 TARAGA, KASUGAI-CITY, AICHI

486-0808 JAPAN

TEL: +81-568-81-7954

[www.tokaimedpro.co.jp/en/](http://www.tokaimedpro.co.jp/en/)

Not available for sale outside Japan  
CE Marking Pending